

LI -BODS-10 : INCUBATOR SHAKER



SALIENT FEATURES :

- Stainless steel inner chamber;
- Four square semicircle transition, convenient to clean;
- Micro-computer control system; LCD display with functions of time setting, parameter memorizing and power reset recovering. (Optional)
- CFC-free refrigeration and gentle breeze circulation.
- Double door seal, tempered glass, silicon rubber seal,
- Inner & Outer door with magnetic rubber seal.
- Inner chamber with light for convenient observation.
- Feature to reduce evaporation/to prevent desiccation of small volume of microbial cultures.

SPECIFICATIONS :

MAKE	LABARD
Model	LI-BODS-10
Temperature	Inner tank made of 304 grade of stainless steel
Shaking speed	Minimum- 30 rpm and Maximum-400 rpm
Accuracy	± 1 rpm (Microprocessor Controlled)
Orbit Diameter	Adjustable
Platform Size	42 x 42 cm
Drive Mechanism	Magnetic drive
Temperature controller	Microprocessor controlled
Lamp	Programmable UV Lamp
Display	LED Display for all parameters
Clamps	100ml,250ml,500ml and two 1000ml
Capacity	4Litre or more capacity flask.

Registered Office:

P-132, kalindi housing estate,
 Ground floor,
 Kolkata-700089.

Phone No.: +91 33 6444 4699/6050 0008

Mobile: +91 9836748235

E-Mail: sales@rdengineering.co.in/labardinstruments@gmail.com

Web: www.labardinstruchem.com www.labardinstruments.co.in www.rdengineering.co.in

Other City Branch: Mumbai